



Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

1 Identification

- **Product identifier**
- **Trade name:** DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab
- **Article number:**
R1100-50, R1100-250, R1100-8-S, R1200-25, R1200-125, R1150, R1101, R1102, R1103, R1104, R1105, R1106, R1107, R1108, R1109
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



Acute Tox. 4 H302 Harmful if swallowed.
 Skin Irrit. 2 H315 Causes skin irritation.
 Eye Irrit. 2B H320 Causes eye irritation.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard statements**
Harmful if swallowed.
Causes skin and eye irritation.
- **Precautionary statements**
Wear protective gloves.
Wash thoroughly after handling.
Do not eat, drink or smoke when using this product.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.

(Contd. on page 2)



Safety Data Sheet acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 1)

Specific treatment (see on this label).

IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

Rinse mouth.

IF ON SKIN: Wash with plenty of water.

Take off contaminated clothing and wash it before reuse.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



Health = 2

Fire = 0

Reactivity = 0

· **HMIS-ratings (scale 0 - 4)**



Health = 2

Fire = 0

Reactivity = 0

· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** This product is a proprietary solution.

· **Dangerous components:** Void

4 First-aid measures

· **Description of first aid measures**

· **After inhalation:** Supply fresh air; consult a physician in case of complaints.

· **After skin contact:** Immediately wash with water and soap and rinse thoroughly.

· **After eye contact:**

Rinsed opened eye for several minutes under running water. If symptoms persist, consult a doctor.

· **After swallowing:** Do not induce vomiting; immediately call for medical help.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 2)

- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
This product is not flammable.
Use fire fighting measures that suit the environment.
No unusual fire or explosion hazards noted.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Dispose contaminated material as waste according to item 13.
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** No special precautions are necessary if used correctly.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special measures required.
- **Information about storage in one common storage facility:** No special measures required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 3)

- Control parameters
Components with limit values that require monitoring at the workplace:
Additional information:
Exposure controls
Personal protective equipment:
General protective and hygienic measures:
Protection of hands:
Penetration time of glove material:
Eye protection:
Body protection:

9 Physical and chemical properties

- Information on basic physical and chemical properties
General Information
Appearance:
Form: Liquid
Color: Colorless
Odor: No information available
Odor threshold: Not determined.
pH-value: 5.0-7.0
Change in condition
Melting point/Melting range: Undetermined.
Boiling point/Boiling range: Undetermined.
Flash point: Not applicable.
Flammability (solid, gaseous): Not applicable.
Ignition temperature:
Decomposition temperature: Not determined.
Auto igniting: Product is not selfigniting.
Danger of explosion: Product does not present an explosion hazard.
Explosion limits:
Lower: Not Applicable

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 4)

Table with 2 columns: Property and Value. Rows include Upper, Vapor pressure, Density, Relative density, Vapor density, Evaporation rate, Solubility in / Miscibility with Water, Partition coefficient (n-octanol/water), Viscosity (Dynamic, Kinematic), Solvent content (Organic solvents, Other information).

10 Stability and reactivity

- Reactivity Strong oxidizing agents and acids.
Possibility of hazardous reactions Strong oxidizing agents and acids.
Conditions to avoid Strong oxidizing agents and acids.
Incompatible materials: Strong oxidizing agents and acids.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:
Primary irritant effect:
on the skin: Irritant to skin and mucous membranes.
on the eye: Irritating to eyes.
Sensitization: No sensitizing effects known.
Additional toxicological information:

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 5)

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
· Aquatic toxicity: No further relevant information available.
· Persistence and degradability No further relevant information available.
· Behavior in environmental systems:
· Bioaccumulative potential No further relevant information available.
· Mobility in soil No further relevant information available.
· Results of PBT and vPvB assessment
· PBT: Not applicable.
· vPvB: Not applicable.
· Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
· Recommendation: Do not dispose with household garbage.
· Uncleaned packagings:
· Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.

14 Transport information

Table with 2 columns: Hazard/Regulation and Status. Rows include UN-Number, DOT, ADN, IMDG, IATA; UN proper shipping name; Transport hazard class(es); Packing group; and Environmental hazards.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 6)

Special precautions for user Not applicable.

Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code Not applicable.

UN Model Regulation: not regulated

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Section 355 (extremely hazardous substances):

None of the ingredients is listed.

Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

TSCA (Toxic Substances Control Act):

None of the ingredients is listed.

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS07

Signal word Warning

Hazard statements

Harmful if swallowed.

(Contd. on page 8)



Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 7)

Causes skin and eye irritation.

· **Precautionary statements**

Wear protective gloves.

Wash thoroughly after handling.

Do not eat, drink or smoke when using this product.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Specific treatment (see on this label).

IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

Rinse mouth.

IF ON SKIN: Wash with plenty of water.

Take off contaminated clothing and wash it before reuse.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 01/16/2017 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

(Contd. on page 9)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 9/9

Safety Data Sheet
acc. to OSHA HCS

Printing date 01/16/2017

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab

(Contd. of page 8)

PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit
Acute Tox. 4: Acute toxicity – Category 4
Skin Irrit. 2: Skin corrosion/irritation – Category 2
Eye Irrit. 2B: Serious eye damage/eye irritation – Category 2B

US